

# BI+ ECO BS300 Vein Finder




## Parameter

- Light Type: Near-infrared Light
- Best Projection Distance: 200±20cm
- Dimension: 24\*60\*60cm(L\*W\*H)
- Battery Working Time: 4 hours
- Operating Mode: Handheld or suspension

## Visualization Modes

- **Brightness**  
Adjustable brightness allows clinicians customizing the image brightness base on the environment light.
- **Colors**  
Switchable colors can fit different needs from the nurses.
- **Inverse**  
Clinicians can highlight the veins or highlight the surrounding tissue.

## Optional Stand



## Application Images



## Application Department



**Wuxi Biomedical Technology Co., Ltd.**

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cncbiomedical.com

Mobile phone: 86-13382885910

[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.